



## Interprosthetic and Mono RescueSleeve® Custom-Made Prostheses

Presented by:



**Waldemar Link GmbH & Co. KG**

Barkhausenweg 10 · 22339 Hamburg, Germany  
Tel.: +49 40 53995-0 · Fax: +49 40 5386929  
info@linkhh.de · [www.linkorthopaedics.com](http://www.linkorthopaedics.com)

# Interprosthetic and Mono RescueSleeve<sup>®</sup> Custom-Made Prostheses

## **04 Description**

Features of the Interprosthetic Twin Sleeves

## **05 Basic Examples**

Twin RescueSleeve<sup>®</sup> with Modular Bridge Coupling for interprosthetic fractures

One-sided RescueSleeve<sup>®</sup> as Custom (A) or MEGASYSTEM-C<sup>®</sup> (B) Coupling

One-sided RescueSleeve<sup>®</sup> with Stem or Hinge

## **06 Surgical Technique**

## **08 Case Examples – Lower Extremities**

## **22 Case Examples – Upper Extremities**

## **24 Additional Instruments**

## **25 Literature**



### Features of the Interprosthetic Twin RescueSleeve®

- Sleeves to take the proximal and distal stem
- Strength proven <sup>6,7</sup> modular bridge
- Stem fixation screws
- Modular bridge „Pocket Coupling“ and security screws

The increasing use of stem prostheses in arthroplasty means that stems of hip and/or knee implants quite often have to be coupled in the femur. According to Soenen<sup>1</sup>, there is a dramatic increase in the risk of interprosthetic fracture, when the two stem tips are less than 100 mm apart. This scenario was tested in a four point bending test.

Weiser et al<sup>2</sup> conducted similar tests and concluded that the distance between the stem tips has scarcely any influence on the risk of fracture, whereas the bone quality of the cortex is the decisive factor.

No matter what the main factor is, good permanent stabilization of such interprosthetic fractures can be achieved with LINK® RescueSleeve®.

Sleeve couplings by their nature do not require narrow tolerances and they have proved very successful for stem/stem coupling<sup>3</sup>. This applies both to stems of LINK implants and to combinations of LINK and other makes of prosthesis (stem dia. must be available).

Essentially, these sleeves can be either single-sided or double-sided (twin sleeve). In the case of the single-sided sleeve, the component that is in contact with the sleeve can comprise an intramedullary stem or a joint component that anchors the prosthesis, which following coupling is held by the sleeve in the medullary canal, or connects it with it opposite side modularly to another joint prosthesis. The double-sided sleeve connects two opposite prostheses stems usually following an interprosthetic fracture, often after a failed osteosynthesis.

The RescueSleeve® can be designed so that an anatomical angle (varus/valgus) is provided in the coupling of the connecting components.

The sleeves are manufactured, for example, in the following sizes:

50 mm depths	8 fixation screws
60 mm depths	12 fixation screws
80 mm depths	16 fixation screws

These sleeve connections always entail a certain bone loss of 130 - 170 mm, but the coupling is so strong that the held prosthesis stem usually would bend or fracture before the connection fails.

The pocket-coupling between the sleeves has a length of 50 mm and is usually applicable.

In few cases of short interprosthetic space the CS (Cardan-Short) coupling has to be employed as that has only 30 mm of length.


Most importantly, the joint regions in the knee or hip of the prosthesis in situ are not impaired by the intervention. Tests on the stability of the sleeve coupling were performed by Professor Morlock at Hamburg University of Technology (TUHH).<sup>4,6</sup>

To create the sleeve/stem coupling, the sleeve is first filled with bone cement before the male in situ component (stem) is pushed into the sleeve while the cement is still soft. Primary fixation of the male component is then achieved with the circumferentially arranged fixation screws. Once the bone cement has hardened, the result is a stable, loadable connection between the in situ prosthesis stem and the attached sleeve (see surgical technique).


In case of the twin RescueSleeve®, the opposing fracture segment is attached to the adjacent sleeve in the same way before both sleeve components are inter-connected via the 'Pocket Coupling' and secured with two assembly screws.

**These implants are only available as custom-made products.**

**Twin RescueSleeve® with modular bridge coupling for interprosthetic fractures**




**Mono RescueSleeve® as custom (A) or MEGASYSTEM-C® (B) coupling**




**A**

**B**

**Mono RescueSleeve® with stem or hinge-connection**




**Mono RescueSleeve® with stem or hinge-connection**



**A**


**B**


**Mono RescueSleeve® attached to Endo-Model® Knee prosthesis with mounting screws to secure a circular stem**



We suggest our drill sleeve (15-8006/03) and our screwdriver (10-5373/02).

- for drilling the tibia component in situ
- to achieve sufficient fixation of the screws





Following site approach the resection length is determined and resection markings are applied to the bone using the coupled sleeve for templating.



Then the proximal and distal stem portion of the prostheses in situ are freed of bone under observation of the resection marks.

Temporary protective cerclages or bone clamps must be used during this procedure.

(see pages 07 and 20)



Next step is a trial positioning. First the stem fixation screws are introduced into the respective bores but only as far to allow the stem to be introduced. It has to be observed that the screws of proximal and distal sleeve might be of different length, due to the variations of proximal and distal stem diameters.


Consequently the screws must not be mixed and only used for their destined sleeve. After provisional fixation of a stem at one side the opposite stem is as well provisionally fixed in the opposite sleeve under observation of the correct rotational position.



Thereafter the rotational position is marked with a sterile pen at the stem and correspondingly the edge of the opening of the second sleeve.




Leg length and rotation are controlled with both sleeves provisionally connected at their coupling ends.



The stem fixation screws to be actively used are left in place (a minimum of at least one pair screws near to the coupling and one pair at the sleeve opening is necessary).



In each of the sleeves two screws holes next to the coupling should be left blank to allow the air to exhaust, when the cement is introduced, preferably performed with a cement gun having a nozzle.



Then the second stem is introduced in the cement filled sleeve of the second component, under observation of the rotational marks and the screws are alternately tightened as well before cement hardens.

Following hardening of the cement both Rescue Sleeves are coupled in the pocket coupling and secured with the two assembly screws having their UHMWPE lock bolts in place in the screws threaded portion.

**Note:** The fixation screws have a 2,5 mm hex head.

Finally the stem of the first component is pushed the whole way into the cement-filled sleeve and the fix screws are alternately tightened while the cement is still soft. Consecutively hardening of the bone cement has to be waited for.

## Titanium Cerclage Band with Lock






## Case 1


### Interprosthetic fracture due to patient downfall - Twin RescueSleeve®.




Interprosthetic fracture with pseudoarthrosis and fracture of osteosynthesis plate.



Technical proposal



„Pocket Coupling“ disassembled



The implant




Post-OP X-ray




## Case 2

## Interprosthetic fracture due to patient downfall - Twin RescueSleeve®.




## Interprosthetic fracture between hip and knee prostheses



## Technical proposal



### The implant




Post-OP X-ray

## Case 3


Interprosthetic fracture due to patient downfall - Twin RescueSleeve® and coupling to MEGASYSTEM-C® components.




Interprosthetic fracture




Secondary dislocation



Technical proposal



Intraoperative



Intraoperative X-ray



One week Post-OP

## Case 4


### Interprosthetic fracture due to patient downfall - Twin RescueSleeve®.



Interprosthetic fracture



Pre-OP Planning



Technical proposal




Intraoperative



Intraoperative X-ray



The implant



Post-OP X-rays





## Case 6


RescueSleeve® after recurrent periprosthetic fractures. One-sided RescueSleeve® to Wagner stem with modular connection to MEGASYSTEM-C® Rotational Knee with modular proximal tibia replacement.



Proximal:  
Wagner revision  
stem



Distal: Knee TKA  
Recurrent peri-  
prosthetic fractures



## Technical proposal




Distal: Knee TKA  
Recurrent peri-  
prosthetic fractures



### The implant



### Post-OP X-rays



## Case 7


One-sided modular RescueSleeve® with hip-joint replacement. PowerLock MP®-coupling allows for 5° increments of antetorsion adjustment.



Pre-OP X-rays



Pre-OP Planning

**LINK**  **Sonderanfertigung / custom made prosthesis**  
gemäß HL 9343/00M, Änderung VIII 1. aus: HL 9343/00M, Januar VIII

Station / date \_\_\_\_\_

Operator / operating surgeon \_\_\_\_\_

Patient / patient \_\_\_\_\_

Kopfteil		Größe
Art.-Nr.	Umsatz	
1101-01	1101-01	1101-01
1101-02	1101-02	1101-02
1101-03	1101-03	1101-03
1101-04	1101-04	1101-04
1101-05	1101-05	1101-05
1101-06	1101-06	1101-06
1101-07	1101-07	1101-07
1101-08	1101-08	1101-08
1101-09	1101-09	1101-09
1101-10	1101-10	1101-10
1101-11	1101-11	1101-11
1101-12	1101-12	1101-12
1101-13	1101-13	1101-13
1101-14	1101-14	1101-14
1101-15	1101-15	1101-15
1101-16	1101-16	1101-16
1101-17	1101-17	1101-17
1101-18	1101-18	1101-18
1101-19	1101-19	1101-19
1101-20	1101-20	1101-20
1101-21	1101-21	1101-21
1101-22	1101-22	1101-22
1101-23	1101-23	1101-23
1101-24	1101-24	1101-24
1101-25	1101-25	1101-25
1101-26	1101-26	1101-26
1101-27	1101-27	1101-27
1101-28	1101-28	1101-28
1101-29	1101-29	1101-29
1101-30	1101-30	1101-30
1101-31	1101-31	1101-31
1101-32	1101-32	1101-32
1101-33	1101-33	1101-33
1101-34	1101-34	1101-34
1101-35	1101-35	1101-35
1101-36	1101-36	1101-36
1101-37	1101-37	1101-37
1101-38	1101-38	1101-38
1101-39	1101-39	1101-39
1101-40	1101-40	1101-40
1101-41	1101-41	1101-41
1101-42	1101-42	1101-42
1101-43	1101-43	1101-43
1101-44	1101-44	1101-44
1101-45	1101-45	1101-45
1101-46	1101-46	1101-46
1101-47	1101-47	1101-47
1101-48	1101-48	1101-48
1101-49	1101-49	1101-49
1101-50	1101-50	1101-50
1101-51	1101-51	1101-51
1101-52	1101-52	1101-52
1101-53	1101-53	1101-53
1101-54	1101-54	1101-54
1101-55	1101-55	1101-55
1101-56	1101-56	1101-56
1101-57	1101-57	1101-57
1101-58	1101-58	1101-58
1101-59	1101-59	1101-59
1101-60	1101-60	1101-60
1101-61	1101-61	1101-61
1101-62	1101-62	1101-62
1101-63	1101-63	1101-63
1101-64	1101-64	1101-64
1101-65	1101-65	1101-65
1101-66	1101-66	1101-66
1101-67	1101-67	1101-67
1101-68	1101-68	1101-68
1101-69	1101-69	1101-69
1101-70	1101-70	1101-70
1101-71	1101-71	1101-71
1101-72	1101-72	1101-72
1101-73	1101-73	1101-73
1101-74	1101-74	1101-74
1101-75	1101-75	1101-75
1101-76	1101-76	1101-76
1101-77	1101-77	1101-77
1101-78	1101-78	1101-78
1101-79	1101-79	1101-79
1101-80	1101-80	1101-80
1101-81	1101-81	1101-81
1101-82	1101-82	1101-82
1101-83	1101-83	1101-83
1101-84	1101-84	1101-84
1101-85	1101-85	1101-85
1101-86	1101-86	1101-86
1101-87	1101-87	1101-87
1101-88	1101-88	1101-88
1101-89	1101-89	1101-89
1101-90	1101-90	1101-90
1101-91	1101-91	1101-91
1101-92	1101-92	1101-92
1101-93	1101-93	1101-93
1101-94	1101-94	1101-94
1101-95	1101-95	1101-95
1101-96	1101-96	1101-96
1101-97	1101-97	1101-97
1101-98	1101-98	1101-98
1101-99	1101-99	1101-99
1101-100	1101-100	1101-100

Waldemar Link GmbH & Co. KG  
Implantate und chirurgische Instrumente  
Birkhäuserweg 12 · D-22329 Hamburg

- unmissverständlich  
an alle  
- alle Maße in mm  
all measurements in mm

Für die Lipung der Sonderanfertigung ist alleine der bestellende Arzt verantwortlich.  
The ordering physician is solely responsible for determining the suitability of this custom-made product.

Technical proposal



The implant

left:  
The sleeve with push-through stem  
and temporary closure screw ...



right:  
... attached to head-neck  
segment.




Post-OP X-ray

## Case 8


### One-sided RescueSleeve® with proximal femur and hip replacement.



Pre-OP X-rays



Technical proposal




Pre-OP Planning

Connection to a femoral long stem of an Endo-Model® Rotational Knee Prosthesis in situ with lateral stem migration. Distal plasmoytome with low activity and osteolysis.

The RescueSleeve® with 'Power Lock' connection to a proximal head-neck component of the LINK® MP® Reconstruction Hip System.

The proximal part of the long stemmed knee could be coupled via RescueSleeve® with the MP® head-neck because LINK custom long stems standard-wise receive the 'Power Lock' connection for possible future component attachments.

The implant



Post-OP X-rays






## Case 10


MEGASYSTEM-C® distal femur replacement having one-sided modular RescueSleeve®.  
Hinged version. Tibial monoblock prosthesis, including RescueSleeve® and fixation screws.



Pre-OP X-ray



## Technical proposal



## Pre-OP Planning



## The implant



### Post-OP X-rays












## Case 13

**Twin RescueSleeve® for interprosthetic fracture between the cemented SP II® stem and a cemented knee prosthesis (competitor implant) in following fracture of the bone plate.**



Pre-OP X-rays



Technical proposal



Pre-OP Planning

Trumpet-Type distal sleeve for optimal seating



The implant



Post-OP X-rays



## Case 14

**Twin RescueSleeve®** in case of repeated fracture of the femur after treatment with bone plate after three months.




Pre-OP



Post-OP



Second fracture  
after 3 months



Technical proposal



LINK proposal



Intraoperative



The implant



Post-OP X-ray



Proximal MEGA-C  
implant components






## Case 16


Humeral twin RescueSleeve® at the right side uniting a shoulder and an elbow prosthesis.



Pre-OP X-ray



## Technical proposal




Pre-OP with proposal



Intraoperative




The implant with CS (Cardan Short) coupling



Post-OP X-ray

### Thabe Titanium Cerclage Band



**63-4300/02** Titanium Cerclage Band with Lock



### **72-1077** Swivel-Jaw Bone Holding Clamp




### **10-5373/01** Hex screwdriver



### **16-3290/00** Cross slot screwdriver



### **15-8006/03** Drill sleeve




### **10-5373/02** Hex screwdriver



- 1 - Soenen, Marc et al „Stemmed TKA in a Femur with a Total Hip Arthroplasty. Is there a safe distance between the stem tips?“, Journ. of Arth., 28 (2013) 1437-1445
- 2 - Weiser, L. et al „The role of interprosthetic fractures of the femur“, JBJS Vol 96-B, No 10, Oct 2014, PP 1378-1384
- 3 - Citak, Mustafa et al „Treatment of interprosthetic femoral fractures with an interposition prosthesis“, Acta Orthopeda 2013, 84 (3): 326-327
- 4 - Weiser, L. (UKE), Morlock M. (TUHH) et al „Biomechanische Untersuchung zum Einfluss des Abstandes zwischen zwei intramedullären Implantaten sowie der Knochendichte und Knochendicke bei interprosthetischen Femurfrakturen“, UKE und TUHH Hamburg, 2013
- 5 - Patel, Nirav K. et al „Custom-made Cement-Linked Mega Prostheses: A Salvage Solution for Complex Periprosthetic Femoral Fractures“, The Journ. of Arth. 29 (2014) 204-209
- 6 - Duda, Georg et al „Dynamic Examination of an Arthrodesis Nail's Taper Connection“, TU Hamburg/Harburg, Dept. of Biomechanics, Prof. Dr. E. Schneider, 02/94
- 7 - Weiser, L. et al „Interposition sleeve as treatment option for interprosthetic fractures of the femur; a biomechanical in vitro assessment“, Intern. Orthop (SICOT), DOI 10.1007/s00264-015-2788-5

## Additional Information

Catalogs on request: E-mail [customer@linkhh.de](mailto:customer@linkhh.de)



255\_Cerclagen



746\_Endo-Model® Knee Fusion Nail SK



Waldemar Link GmbH & Co. KG

Barkhausenweg 10 • 22339 Hamburg, Germany

Phone +49 40 53995-0 • [info@linkhh.de](mailto:info@linkhh.de)

[www.linkorthopaedics.com](http://www.linkorthopaedics.com)

