

**deSouther**  
MEDICAL



**osteodrive<sup>®</sup>**

[www.de-soutter.com](http://www.de-soutter.com)



# osteodrive<sup>®</sup>

- The osteodrive™ electric small bone range has been carefully designed to provide seamless system intelligence and customisable instrument performance.
- The osteodrive™ unique attachment recognition simplifies user operation, whilst also improving hand piece balance and size.
- Fully customisable and adjustable performance settings allow the exact requirement of the surgeon to be met.
- Intelligent irrigation recognition means no set up required, and unique magnetic nozzle mounts provide stable yet flexible use.
- Quick and simple lever removal ensures users can quickly assemble the system for any mode of operation.

# MCI MODULAR HANDPIECES

The De Soutter Medical osteodrive™ is a versatile modular electric instrument system for small bone surgery. The handpieces are powered by a high performance, dual socket power supply providing greater surgical control and precision.



MCI-270  
With mini sagittal  
saw attachment



MCI-270  
With rotary bur and  
detachable lever

MCI-270  
With  
wire driver  
attachment




Lever control Order code 18710  
Easily detachable lever  
for MCI-270 handpieces



Lever extends for  
optimal finger  
positioning

Handpiece	Control	Order code	Speed	Power supply
MCI-270	Foot pedal or lever	1292004	0 - 35000 rpm	Via 4 metre cable



## MCI DRIVES AND ATTACHMENTS

	Order code	Model	Description	Speed rpm
<b>DRILL ATTACHMENTS</b>				
	18490	DI-270	0.5 - 4.0 mm (5/32") Keyed	0 - 1450
	19240	DI-270	Inline A.O.	0 - 1450
	19320	DI-270	Offset A.O. 2.2 mm cannulation	0 - 1500
	19580	DI-270	Offset Jacobs (5/32")	0 - 1500
<b>HIGH SPEED BURS</b>				
	18520	BI-270	Accepts all 2.38 mm (3/32") burs Uses De Soutter Medical R Series burs	0 - 35000
	19460	BI-277	Low Speed Bur. Uses R Series burs (BI-2xx supplied with 19470 attachment holder)	0 - 9000
	18720	-	Standard bur guard	-
	18730	-	Long bur guard	-
	18740	-	Extra long bur guard	-
	18590	QI-270	Medium speed drill ISO 1797 Type1	0 - 5000
	18540	BI-278	Metal cutter - includes splash hood Uses De Soutter Medical R60 Series burs	0 - 17500
	18480	EI-270	E-Drive	0 - 35000
	14940	AV-221	200 Contra-angle rotary bur 1:1 Used with E-Drive 18480 Accepts 2.35mm RA shanks	-
	9770	AV-220	Standard length 20 ° rotary bur 1:1	-
	9780	AV-220	Extended length 200 rotary bur 1:1 Used with E-Drive Order code 18480 Accepts 2.35mm (3/32") R Series burs	-
<b>SAGITTAL SAWS</b>				cpm
	18510	SI-270	Mini sagittal saw Uses S86 saw blades	0 - 27000
	18460	NI-270	Large sagittal saw Uses S88 saw blades	0 - 22000
<b>RECIPROCATING SAW</b>				
	18500	CI-270	Small reciprocating saw Uses S83 saw blades	0 - 14500
	19410	CI-271	Small reciprocating saw Uses S92 blades	0 - 14500
<b>OSCILLATING SAW</b>				
	19820	OI-270	Oscillating saw Uses S86 type blades	0 - 22000
	19920	OI-271	Oscillating saw Uses S86-9XXX type blades	0 - 20000
<b>WIRE DRIVERS</b>				
	18560	WI-271	In-line wire driver Capacity 0.6 - 1.6 mm	0 - 1450
	18550	WI-270	Angled wire driver Capacity 0.6 - 1.6 mm	0 - 1450

Ask for Stericut® catalogue for range of sterile packed saw blades, drill bits, rotary burs and K wires.

# MCU MODULAR HANDPIECES



MCU-470  
With quick  
release  
drill attachment

































MCU-470  
With quick  
release  
sagittal saw  
attachment



Handpiece	Order code	Speed	Cannulation	Power supply
MCU-470	1291994	0 - 1350 rpm	4.2 mm	Via 4 metre cable


## MCU DRIVES AND ATTACHMENTS

	Order code	Model	Chuck description	Speed rpm	Cannulation Ø mm
<b>DRILL ATTACHMENTS</b>					
	16670	DU-470	0.5 - 4.0 mm (5/32") Keyed	0 - 1350	4.1
	16680	DU-470	0.5 - 7.4 mm (1/4") Keyed	0 - 1350	4.2
	16690	DU-470	0.3 - 3.0 mm (1/8") Keyless	0 - 1350	3.0
	16700	DU-470	0.3 - 7.4 mm (1/4") Keyless	0 - 1350	4.2
<b>QUICK RELEASE DRILLING ATTACHMENTS</b>					
	16710	DU-470	A.O. Synthes/ASIF	0 - 1350	2.5
	16720	DU-470	Hudson/Zimmer	0 - 1350	4.2
	16730	DU-470	Trinkle	0 - 1350	4.2
	19590	DU-470	DHS/DCS	0 - 1350	4.2
	20260	DU-470	Aesculap	0 - 1350	-
	16740	QU-470	Medium speed drill ISO 1797 Type 1	0 - 4500	-
<b>HIGH SPEED BUR</b>					
	19440	BU-472	High speed bur - Uses R Series burs	0 - 27000	-
			Standard bur guard Order code 18720		
			Long bur guard Order code 18730		
			Extra long bur guard Order code 18740		
<b>REAMING ATTACHMENTS</b>					
	16750	RU-470	0.5 - 7.4 mm (1/4")	0 - 350	4.2
	17750	RU-471	0.5 - 7.4 mm (1/4") - non reversible	0 - 350	4.2
	16760	RU-470	A.O. Synthes	0 - 350	4.2
	17740	RU-471	A.O. Synthes - non reversible	0 - 350	4.2
	16770	RU-470	Hudson/Zimmer	0 - 350	4.2
	19340	RU-470	DHS/DCS	0 - 350	4.2
	17490	RU-470	Hudson	0 - 350	4.2
	17760	RU-471	Hudson - non reversible	0 - 350	4.2
	17500	RU-470	Zimmer	0 - 350	4.2
	17760	RU-471	Zimmer - non reversible	0 - 350	4.2
<b>SAWING ATTACHMENTS</b>					
	18700	SU-471	Button release sagittal saw (non rotating) Uses S86 type blades	0 - 22000	-
	18260	SU-472	Button release sagittal saw (rotating) Uses S86 type blades	0 - 22000	-
	16790	NU-470	Quick release sagittal saw Uses S88-1XX blades	0 - 20500	-
	19520	NU-472	High speed sagittal saw (gold) Uses S88-1XX blades	0 - 24300	-
	18070	NU-471	Quick release sagittal saw Uses S88-1XX & S88-4XX blades	0 - 15000	-
	16800	CU-470	Small reciprocating saw Uses S92 type blades	0 - 8000	-
	16810	CU-471	Large reciprocating saw Uses S22 type blades	0 - 7000	-
<b>PULSE LAVAGE ATTACHMENT</b>					
	16820	LU-470	Pulse lavage attachment Uses P30-7XX series pump kits	0 - 1350 pulses/min	-
<b>TORQUE LIMITERS</b>					
	16860	TU-404	0.4 Nm	0 - 350	-
	16870	TU-408	0.8 Nm	0 - 350	-
	16880	TU-415	1.5 Nm	0 - 350	-
<b>RADIOLUCENT DRILL ADAPTOR</b>					
	16830	TU-470	Accepts A.O. radiolucent drill attachment Order code 17710	0 - 1350	-
<b>WIRE/PIN DRIVERS</b>					
	17180	WU-471	0.7 - 4.0 mm capacity	0 - 1350	4.0
	16840	WU-470	0.6 - 1.8 mm capacity	0 - 1350	1.8
	16850	PU-470	1.8 - 3.2 mm capacity	0 - 1350	3.2

Ask for Stericut® catalogue for range of sterile packed saw blades, drill bits, rotary burs and K wires.

# POWER CONSOLES

The osteodrive power console range provides a step change in performance for the end user




- Simultaneous dual handpiece operation.
- Compatible with a range of systems already on-site, ensuring cost containment and reducing the need for multiple power sources.
- On-site upgrades are easily and quickly performed using the USB port
- Audible safety mode.




### 1 Customisable

200+ user profiles can be stored and are easily accessed  
Multi language options



### 2 Tailored Performance

Adjustable acceleration, braking and torque allows precise user control for every specific case



### 3 Ease of use

Wireless foot pedal control ensures efficient operation  
Large 7" colour touch screen



### 4 Intuitive design

Ports illuminate to clearly indicate tool status to all OR staff  
Attachment recognition communicates tool status



# POWER CONSOLES

## PC-47X POWER CONSOLES



Model	Description	Order code
PC-472	Advanced power console with irrigation	19090



PC-471	Advanced power console	19080
--------	------------------------	-------



PC-470	Basic 'plug & play' power console	18050
--------	-----------------------------------	-------



4 metre cable for connecting handpiece to Power Console  
Order code 18570

Mains power cordsets -3 metres long.  
Please specify -from these options -

	Order code
UK	8290
Euro	8300
US/Canada	8310
Australia/New Zealand	8320
Swiss	14700

## FS-27X FOOT PEDALS

With three options of foot pedal, every function is available from wireless control to full capability across our small bone power range

Description	Model	Order code
Foot pedal - wireless	FS-271	18600
Foot pedal - wired	FS-270	18580
Foot pedal - wired *	FS-272	19910

\* For use with De Soutter Medical hekaDrill range.

Charger for wireless foot pedal  
Order code 19180



FS-271  
Wireless

FS-270  
Wired

## IRRIGATION SYSTEM

The PC-472 provides total versatility for the osteodrive range. A simple front loading irrigation cassette with integrated RFID identification technology guarantees optimal set up on every occasion, avoiding the need for complex set up routines.



MCI-270 attached to PC-472 controller with T20-100 irrigation kit

### T20-100 IRRIGATION KIT



### IRRIGATION NOZZLE




De Soutter Medical's unique magnetic nozzle design provides stable yet flexible use. Magnetised adjustable wings sit firmly to any attachment type ensuring total flexibility and precision when directing fluid. Repositioning of the nozzle is quick and easy due to the unique properties of the magnetic wings.

\*T20-121 For use with De Soutter Medical hekaDrill range.

# STERILISATION TRAYS AND CONTAINERS

## INSERT TRAYS



Set up - optional tray contents  
 1 off MCI-270 handpiece  
 1 off MCU-470 handpiece  
 6 off MCI-270 attachments  
 9 off MCU-470 attachments  
 2 off connection cables

Container and insert tray  
 Order code

18950

18990

Half size wire container  
 and insert tray

Half size  
 Insert tray



1 off MCI-270 handpiece  
 6 off MCI-270 attachments  
 1 off connection cable

19420

19260

## WIRE CONTAINERS

Accepts standard insert tray



Standard wire container  
 Order code 19010

Accepts half size insert tray



Half size wire container  
 Order code 19000

## STERILE FILTERED CONTAINERS

Accepts standard insert tray  
 Standard sterile filtered container  
 Order code 19250

PTFE sterile filtered container  
 Order code 19810

Lid gasket - Order code 670023



PTFE filter cassette  
 Order code 671073



Accepts half size insert tray

Half size Sterile Filtered Container -  
 Order code 19370

Lid gasket - Order code 669903



Filter paper - single use  
 Pack of 100  
 Order code 657243



Security tag  
 Pack of 100  
 Order code 657253



Sterile indicator label  
 Pack of 100 -  
 Order code 657263

**UNITED KINGDOM (INTERNATIONAL HEAD OFFICE)**

De Soutter Medical Limited  
Halton Brook Business Park, Weston Road,  
Aston Clinton, Aylesbury, Bucks HP22 5WF  
Tel: +44 (0) 1296 634 000  
Fax: +44 (0) 1296 634 033  
Email: info@de-soutter.com

**AUSTRALIA**

De Soutter Medical Australia Pty Ltd  
2/12-14 Apollo Drive, Hallam  
Victoria 3803  
Tel: +61 (0) 3 9702 4441  
Fax: +61 (0) 3 9702 4484  
Email: australia@de-soutter.com

**AUSTRIA**

De Soutter Medical Austria  
Zweigniederlassung Österreich  
Dietrichsteingasse 10, A-3400 Klosterneuburg  
Tel: +43 (0) 676 96 71 770  
Fax: +43 (0) 2243 21 656  
Email: austria@de-soutter.com

**BELGIUM**

De Soutter Medical Belgium  
Bessemerstraat 14  
3620 Lanaken  
Tel: +32 (0) 89/47 15 37  
Fax: +32 (0) 89/70 12 19  
Email: belgium@de-soutter.com

**FRANCE**

De Soutter Medical France  
1 252 Avenue du Parc des Expositions  
33260 La Teste de Buch  
Tel: +33 (0) 5 56 54 89 36  
Fax: +33 (0) 9 70 61 37 60  
Email: france@de-soutter.com

**GERMANY**

De Soutter Medical Deutschland  
Bahnhofstraße 4  
D-66625 Nohfelden  
Tel: +49 (0) 68 52-99 12 46  
Fax: +49 (0) 68 52-99 12 47  
Email: deutschland@de-soutter.com

**ITALY**

De Soutter Medical Italia  
Località Fornace SNC  
27022 Casorate Primo - PV  
Tel: +39 (0) 2 9009 4098  
Fax: +39 (0) 2 9009 2673  
Email: italy@de-soutter.com

**NETHERLANDS**

De Soutter Medical Nederland  
Gelderlandhaven 2X  
3433 PG Nieuwegein  
Tel: +31 (0) 85 0491480  
Fax: +31 (0) 85 0491489  
Email: nederland@de-soutter.com

**UNITED STATES OF AMERICA**

De Soutter Medical USA, Inc.  
224 Rolling Hill Road, Suite 12A  
 Mooresville, NC 28117  
Tel: +1 (704) 655 9040  
Fax: +1 (704) 987 2035  
Email: usa@de-soutter.com

**deSoutter**  
MEDICAL

For further information on De Soutter Medical powered instruments and cast saw systems please contact the address in your area.