

## REVISION SOLUTIONS

Multiple configurations  
for hip, knee, extremities  
and spinal procedures





R60-504 - Round head Ø5 mm x 24 mm long  
 R60-506 - Round head Ø7 mm x 27 mm long  
 R60-509 - Round head Ø10 mm x 29 mm long  
 R60-516 - Flame head Ø7 mm x 39 mm long



**HANDPIECE SPECIFICATION**

Model	Order code	Maximum reaming torque	Speed rpm
MBQ-707	1289334	16 Nm	0 - 1350
MBQ-708	1289444	16 Nm	0 - 1350

**ROTARY REVISION BURS**

Model	Order code	Bur length	Speed rpm
BQ-707 Standard	19890	72 mm	0 - 20000
BQ-707 Long	19880	197 mm	0 - 20000
BQ-707 Extra long	19870	273 mm	0 - 20000



Uses R60-5XX Series revision burs

# osteodrive™



MCI fitted with BI-270 attachment



BI-278 metal cutter attachment

## HANDPIECE SPECIFICATION

Model	Order code	Speed rpm
MCI-270	1292004	0 - 35000

## ATTACHMENTS

Model	Order code	Description	Speed rpm
BI-278	18540	Metal cutter. Includes splash hood. Uses De Soutter Medical R60 Series burs	0 - 17500
BI-277	19460	Low speed. Accepts all 2.38 mm (3/32") R Series burs (BI-277 is supplied with Order code 19470 attachment holder)	0 - 9000
BI-270	18520	Accepts all 2.38 mm (3/32") R Series burs.	0 - 35000
Standard bur guard - Order code 18720		Long bur guard - Order code 18730	Extra long bur guard Order code 18740

# hekaDrill



## HANDPIECE SPECIFICATION

Model	Order code	Description	Speed rpm
HDS-121	103347	Electric	200 - 85000



## ATTACHMENTS

Type	70 mm long		80 mm long		100 mm long		150 mm long	
	Model	Order code	Model	Order code	Model	Order code	Model	Order code
Straight Angled	S-C07	103360	S-C08	103371	SC-10	103383	SC-15	103394
	A-C07	103446			AC-10	103451	AC-15	103467
Metal cutter	M-D06	104137						

Please ask for catalogue No. 212 for full details of the De Soutter Medical osteodrive™ and No. 213 for hekaDrill systems.

UNITED KINGDOM (INTERNATIONAL HEAD OFFICE)

De Soutter Medical Limited  
Halton Brook Business Park, Weston Road  
Aston Clinton, Aylesbury, Bucks HP22 5WF  
Tel: +44 (0) 1296 634 000  
Fax: +44 (0) 1296 634 033  
Email: info@de-soutter.com

AUSTRALIA

De Soutter Medical Australia Pty Ltd  
2/12-14 Apollo Drive, Hallam  
Victoria 3803  
Tel: +61 (0) 3 9702 4441  
Fax: +61 (0) 3 9702 4484  
Email: australia@de-soutter.com

AUSTRIA

De Soutter Medical Austria  
Zweigniederlassung Österreich  
Dietrichsteingasse 10, A-3400 Klosterneuburg  
Tel: +43 (0) 676 96 71 770  
Fax: +43 (0) 2243 21 656  
Email: austria@de-soutter.com

BELGIUM

De Soutter Medical Belgium  
Bessemmerstraat 14  
3620 Lanaken  
Tel: +32 (0) 89/47 15 37  
Fax: +32 (0) 89/70 12 19  
Email: belgium@de-soutter.com

FRANCE

De Soutter Medical France  
1 252 Avenue du Parc des Expositions  
33260 La Teste de Buch  
Tel: +33 (0) 5 56 54 89 36  
Fax: +33 (0) 9 70 61 37 60  
Email: france@de-soutter.com

GERMANY

De Soutter Medical Deutschland  
Bahnhofstraße 4  
D-66625 Nohfelden  
Tel: +49 (0) 68 52-99 12 46  
Fax: +49 (0) 68 52-99 12 47  
Email: deutschland@de-soutter.com

ITALY

De Soutter Medical Italia  
Località Fornace SNC  
27022 Casorate Primo - PV  
Tel: +39 (0) 2 9009 4098  
Fax: +39 (0) 2 9009 2673  
Email: italy@de-soutter.com

NETHERLANDS

De Soutter Medical Nederland  
Gelderlandhaven 2X  
3433 PG Nieuwegein  
Tel: +31 (0) 85 0491480  
Fax: +31 (0) 85 0491489  
Email: nederland@de-soutter.com

UNITED STATES OF AMERICA

De Soutter Medical USA, Inc.  
224 Rolling Hill Road, Suite 12A  
 Mooresville, NC 28117  
Tel: +1 (704) 655 9040  
fax: +1 (704) 987 2035  
Email: usa@de-soutter.com

**deSoutter**  
MEDICAL

For further information on De Soutter Medical powered instruments and cast saw systems please contact the address in your area.