

eco  
PULSE

“ THE NEXT GENERATION  
IN PULSE LAVAGE SUSTAINABILITY



[www.de-soutter.com](http://www.de-soutter.com)

**deSoutter**  
MEDICAL



## ECO PULSE™ PULSE LAVAGE SYSTEM



### WHAT SAVINGS DOES THE NEW ECO PULSE OFFER?

#### ECO-FRIENDLY



The ECO Pulse doesn't require batteries, motors or electronics and therefore typically **saves up to 55% in clinical waste**, per procedure.

#### COST



Due to its novel and simple design, the ECO Pulse reduces hospital spend whilst also helping the environment.

#### SPACE



A smaller footprint allows less packaging to be used, this results in a **reduction in valuable storage space by 2.5 times**.

#### USER FRIENDLY

P30 ECO Pulse pump kits simply "click" into any De Soutter Medical rotary handpiece for high performance, instant use.

Our Stryker® compatible P31 series pump kits operate directly with any Stryker® large rotary bone handpiece (System 7 and 8)



P30 Series - secure nozzle connection prevents leaks and breakages


#### PERFORMANCE

Flow rate Variable 0 - 1.0 litres/min

Pulse rate Variable 0 - 1350 pulses/min

\* Based on AB-703 Performance figures dependent on battery type.

# PUMP KITS & NOZZLES



## PUMP KITS

De Soutter Medical Order code	Stryker® Compatible Order code	Type	Supplied nozzle
P30-700	P31-700	General purpose with suction	P30-121
P30-701	P31-701	Hip arthroplasty with suction	P30-121/P30-123
P30-702	P31-702	Hip arthroplasty brush with suction	P30-124/P30-121
P30-704	P31-704	Hip arthroplasty without suction	P30-120/P30-122
P30-705	P31-705	Bone cleaning (high pressure) with suction	P30-125
P30-707	P31-707	Hip arthroplasty (high pressure) with suction	P30-125/P30-123
P30-708	P31-708	Hip arthroplasty (high pressure) without suction	P30-125/P30-122
P30-709	P31-709	Bone cleaning (fan spray) with suction	P30-127
P30-710	P31-710	Hip arthroplasty (fan spray) with suction	P30-123/P30-127

## NOZZLES

Order code	Useable length	Type	
P30-125	6cm	Bone cleaning - high pressure with suction	
P30-127	6cm	Bone cleaning - fan spray with suction	
P30-121	6cm	Shower spray with suction	
P30-120	10cm	Soft tissue without suction	
P30-123	25cm	Femoral canal with suction	
P30-124	25cm	Femoral brush with suction	
P30-122	25cm	Femoral canal without suction	

Suction tubing is supplied with all pump kit types. All sterile packed pump kits and individual nozzles are supplied 5 per box.



**UNITED KINGDOM (INTERNATIONAL HEAD OFFICE)**

De Soutter Medical Limited  
Halton Brook Business Park, Weston Road Aston Clinton,  
Aylesbury, Bucks HP22 5WF  
Tel: +44 (0) 1296 634 000  
Fax: +44 (0) 1296 634 033  
Email: info@de-soutter.com

**AUSTRALIA**

De Soutter Medical Australia Pty Ltd  
2/12-14 Apollo Drive, Hallam  
Victoria 3803  
Tel: +61 (0) 3 9702 4441  
Fax: +61 (0) 3 9702 4484  
Email: australia@de-soutter.com

**AUSTRIA**

De Soutter Medical Austria  
Zweigniederlassung Österreich  
Dietrichsteingasse 10, A-3400 Klosterneuburg  
Tel: +43 (0) 676 96 71 770  
Fax: +43 (0) 2243 21 656  
Email: austria@de-soutter.com

**BELGIUM**

De Soutter Medical Belgium  
Bessemmerstraat 14  
3620 Lanaken  
Tel: +32 (0) 89/47 15 37  
Fax: +32 (0) 89/70 12 19  
Email: belgium@de-soutter.com

**FRANCE**

De Soutter Medical France  
1 252 Avenue du Parc des Expositions  
33260 La Teste de Buch  
Tel: +33 (0) 5 56 54 89 36  
Fax: +33 (0) 9 70 61 37 60  
Email: france@de-soutter.com

**GERMANY**

De Soutter Medical Deutschland Bahnhofstraße 4  
D-66625 Nohfelden  
Tel: +49 (0) 68 52-99 12 46  
Fax: +49 (0) 68 52-99 12 47  
Email: deutschland@de-soutter.com

**ITALY**

De Soutter Medical Italia  
Località Fornace SNC  
27022 Casorate Primo - PV  
Tel: +39 (0) 2 9009 4098  
Fax: +39 (0) 2 9009 2673  
Email: italy@de-soutter.com

**NETHERLANDS**

De Soutter Medical Nederland Gelderlandhaven 2X  
3433 PG Nieuwegein  
Tel: +31 (0) 85 0491480  
Fax: +31 (0) 85 0491489  
Email: nederland@de-soutter.com

**UNITED STATES OF AMERICA**

De Soutter Medical USA, Inc.  
224 Rolling Hill Road, Suite 12A Mooresville, NC 28117  
Tel: +1 (704) 655 9040  
Fax: +1 (704) 987 2035  
Email: usa@de-soutter.com

For further information on De Soutter Medical powered instruments please contact the address in your area

*Stryker® is a registered trademark of Stryker Corporation.*

**deSoutter**  
MEDICAL