



# hekaDrill




## ULTIMATE PRECISION

*The hekaDrill has been developed with innovation and precision as a guiding principle.*

Key capabilities include:

- High speed motor control - up to 85,000 rpm. Allows high torque/low speed control. (min. 200 rpm)
- Advanced motor control ensures comfort for the surgeon during extended use with minimal heat gain.
- Adjustable bur exposure lengths.
- Wireless foot pedal operates using a secure encrypted network, ensuring greater user protection.
- Adjustable acceleration and braking functionality.




*A single rotary bur allows 5 possible exposure positions offering total flexibility and reduced inventory costs.*

## ELECTRIC HANDPIECE AND ATTACHMENTS

HIGHSPEED HANDPIECE		Order code	Model	Description	Speed rpm
	hekaDrill	103347	HDS-121	Cable length 4m	200 - 85,000
	Lever Control	103352	HDS-122	Lever Control	200 - 85,000
	Lever	104184	Lever		
ATTACHMENTS STRAIGHT		Order code	Model	Length - mm	Bur size - mm
	S-C07	103360	S-C07	70	2.38
	S-C08	103371	S-C08	80	2.38
	S-C10	103383	S-C10	100	2.38
	S-C15	103394	S-C15	150	2.38
	S-D09	103409	S-D09	90	3.18
	S-D14	103414	S-D14	140	3.18
	S-D26	103438	S-D26	260	3.18
SLIMLINE		Order code	Model	Length - mm	Bur size - mm
	S-C07S	104035	S-C07S	70	2.38
	S-C08S	104043	S-C08S	80	2.38
	S-C10S	104056	S-C10S	100	2.38
	S-C12S	104312	S-C12S	120	2.38
	S-C15S	104062	S-C15S	150	2.38
ANGLED		Order code	Model	Length - mm	Bur size - mm
	A-C07	103446	A-C07	70	2.38
	A-C10	103451	A-C10	100	2.38
	A-C15	103467	A-C15	150	2.38
	A-D14	103472	A-D14	140	3.18
SLIMLINE		Order code	Model	Length - mm	Bur size - mm
	A-C07S	104070	A-C07S	70	2.38
	A-C10S	104081	A-C10S	100	2.38
	A-C12S	104320	A-C12S	120	2.38
	A-C15S	104097	A-C15S	150	2.38
ANGLED MODULAR		Order code	Model	Description	
	A-M01	104159	A-M01	Angled modular attachment	
CRANIOTOME		Order code	Model	Length - mm	Bur size - mm
	C-C09	103508	C-C09	90 (Small)	2.38
	C-C10	103480	C-C10	100 (Medium)	2.38
	C-D12	103495	C-D12	120 (Large)	3.18
METAL CUTTER		Order code	Model	Speed rpm	Chuck type
	M-DO6	104137	M-DO6	-	3.18
PERFORATOR ATTACHMENT		Order code	Model	Speed rpm	Chuck type
	P-P01	103519	P-P01	1,000	Hudson
	Used with HDS-121 at 72,000 rpm				
	Special order only				
	Locking spanner for HDS-1xx supplied with handpiece	103686			
	Ask for catalogue 2677 for our range of rotary burs				



Description	Model	Order code
Console - less irrigation	PC-120	103305
Console - with irrigation	PC-121	103318



The advanced hekaDrill power console is available with or without irrigation and offers numerous features:

- Simultaneous dual handpiece operation.
- Adjustable acceleration and braking functionality.
- Integrated and customisable irrigation control with anti drip mode.
- 200+ user profiles can be stored and easily accessed.
- Rings illuminate to indicate port status for ease of use.

# IRRIGATION AND FOOT PEDALS

## STERILE IRRIGATION KIT



Sterile irrigation kit.  
Box of 5.  
Order code: S101461

## FOOT PEDALS



- Both wired and wireless foot pedal available.
- Wireless foot pedal operates using a secure encrypted network for complete protection.
- A or B tool selection.
- Foot pedal underlighting to allow for easy detection.



4 metre cable for connecting and charging foot pedals via power console.

Order code: 103597

	Model	Order code
Wired	MFP-120	103556
Wireless	MFP-121	103562

## POWER SUPPLY

Mains power cable 3 metres long. Please specify from below.

UK	Order code: 8290
Euro	Order code: 8300
USA/Canada	Order code: 8310
Australia/New Zealand	Order code: 8320
Japan	Order code: 8330
Swiss	Order code: 14700



HANDPIECE	Order code	Model	Speed rpm
-----------	------------	-------	-----------



103339	HDS-111	100,000
--------	---------	---------

The HDS-111 handpiece accepts all straight and angled attachments and craniotome

PERFORATOR	Order code	Model	Speed rpm	Operating pressure
------------	------------	-------	-----------	--------------------



103524	P-P02	850	7 bar
--------	-------	-----	-------

Uses Hudson chuck type

FOOT CONTROL	Order code	Model	Foot pedal inlet connection	Foot pedal outlet connection
--------------	------------	-------	-----------------------------	------------------------------



14900	FV-220	Integral 3 metre hose with MA7 probe	Accepts 'wash assist' hose from sizes below
-------	--------	--------------------------------------	---



3 metre length - Order code 632333  
 4 metre length - Order code 632343  
 5 metre length - Order code 632353



Wash cap for use with HDS-111/ P-P02 Order code 9950

# STERILISATION TRAYS AND CONTAINER

## STERILISATION TRAYS AND CONTAINER



### ELECTRIC INSTRUMENT SYSTEM

Order code: 103610  
Includes lid, container, electric tray and attachment tray.

Handpiece tray  
Optional tray contents -  
• 1 off HDS-12X handpiece  
• 1 off P-PO1 perforator  
• 1 off electric cable  
• 1 off spanner

Attachment tray  
Optional tray contents -  
• 5 off attachments  
• 1 off spanner

### PNEUMATIC INSTRUMENT SYSTEM

Order code 103629  
Includes lid, container, pneumatic tray and attachment tray.

Handpiece tray  
Optional tray contents -  
• 1 off HDS-111 handpiece  
• 1 off P-PO2 perforator  
• 1 off pneumatic hose  
• 1 off spanner

Attachment tray  
Optional tray contents -  
• 5 off attachments  
• 1 off spanner

### CONTAINER AND TRAYS ONLY

Wire container and lid  
Order code: 103944

Tray for electric system  
Order code: 103953

Tray for pneumatic system  
Order code: 103989

Tray for attachments  
Order code: 103966

#### **UNITED KINGDOM (INTERNATIONAL HEAD OFFICE)**

De Soutter Medical Limited  
Halton Brook Business Park, Weston Road,  
Aston Clinton, Aylesbury, Bucks HP22 5WF  
Tel: +44 (0) 1296 634 000  
Email: info@de-soutter.com

#### **AUSTRALIA**

De Soutter Medical Australia Pty Ltd  
2/12-14 Apollo Drive, Hallam  
Victoria 3803  
Tel: +61 (0) 3 9702 4441  
Email: australia@de-soutter.com

#### **AUSTRIA**

De Soutter Medical Austria  
Hauptstrasse 122-124/Büro 01  
3420 Kritzensdorf  
Tel: +43 (0) 2243 239000  
Email: austria@de-soutter.com

#### **BELGIUM**

De Soutter Medical Belgium  
Bessemmerstraat 14  
3620 Lanaken  
Tel: +32 (0) 89/47 15 37  
Email: belgium@de-soutter.com

#### **FRANCE**

De Soutter Medical France  
1 252 Avenue du Parc des Expositions  
33260 La Teste de Buch  
Tel: +33 (0) 5 56 54 89 36  
Email: france@de-soutter.com

#### **GERMANY**

De Soutter Medical Deutschland  
Bahnhofstraße 4  
D-66625 Nohfelden  
Tel: +49 (0) 68 52-99 12 46  
Email: deutschland@de-soutter.com

#### **ITALY**

De Soutter Medical Italia  
Località Fornace SNC  
27022 Casorate Primo - PV  
Tel: +39 (0) 2 9009 4098  
Email: italy@de-soutter.com

#### **NETHERLANDS**

De Soutter Medical Nederland  
Gelderlandhaven 2X  
3433 PG Nieuwegein  
Tel: +31 (0) 85 0491480  
Email: nederland@de-soutter.com

#### **UNITED STATES OF AMERICA**

De Soutter Medical USA, Inc.  
224 Rolling Hill Road, Suite 12A  
 Mooresville, NC 28117  
Tel: +1 (704) 655 9040  
Email: usa@de-soutter.com

**deSoutter**  
MEDICAL

Distributed by De Soutter Medical, manufactured by Zethon