

deSoutter
MEDICAL



hekaDrill


www.de-soutter.com

ULTIMATE PRECISION

The hekaDrill has been developed with innovation and precision as a guiding principle.

Key capabilities include:

- High speed motor control - up to 85,000 rpm. Allows high torque/low speed control. (min. 200 rpm)
- Advanced motor control ensures comfort for the surgeon during extended use with minimal heat gain.
- Adjustable bur exposure lengths.
- Wireless foot pedal operates using a secure encrypted network, ensuring greater user protection.
- Adjustable acceleration and braking functionality.




A single rotary bur allows 5 possible exposure positions offering total flexibility and reduced inventory costs.

ELECTRIC HANDPIECE AND ATTACHMENTS

HIGHSPEED HANDPIECE	Order code	Model	Description	Speed rpm
	103347	HDS-121	Cable length 4m	200 - 85,000
	103352	HDS-122	Lever Control	200 - 85,000
	104184	Lever		
ATTACHMENTS STRAIGHT	Order code	Model	Length - mm	Bur size - mm
	103360	S-C07	70	2.38
	103383	S-C10	100	2.38
	103394	S-C15	150	2.38
	103409	S-D09	90	3.18
	103414	S-D14	140	3.18
	103438	S-D26	260	3.18
SLIMLINE				
	104035	S-C07S	70	2.38
	104056	S-C10S	100	2.38
	104312	S-C12S	120	2.38
	104062	S-C15S	150	2.38
ANGLED				
	103446	A-C07	70	2.38
	103451	A-C10	100	2.38
	103467	A-C15	150	2.38
	103472	A-D14	140	3.18
SLIMLINE				
	104070	A-C07S	70	2.38
	104081	A-C10S	100	2.38
	104320	A-C12S	120	2.38
	104097	A-C15S	150	2.38
ANGLED MODULAR				
	104159	A-M01	Angled modular attachment	
FOOTED ATTACHMENT				
	103508	C-C09	90 (Small)	2.38
	103480	C-C10	100 (Medium)	2.38
	103495	C-D12	120 (Large)	3.18
METAL CUTTER	Order code	Model	Speed rpm	Chuck type
	104137	M-D06	-	3.18
PERFORATOR ATTACHMENT				
	103519	P-P01	Hudson Used with HDS-121 at 72,000 rpm	1,000
			Special order only	
	103686		Locking spanner for HDS-1xx supplied with handpiece	
Ask for catalogue 2677 for our range of rotary burs				



Description	Model	Order code
Console - less irrigation	PC-120	103305
Console - with irrigation	PC-121	103318



The advanced hekaDrill power console is available with or without irrigation and offers numerous features:

- Simultaneous dual handpiece operation.
- Adjustable acceleration and braking functionality.
- Integrated and customisable irrigation control with anti drip mode.
- 200+ user profiles can be stored and easily accessed.
- Rings illuminate to indicate port status for ease of use.

IRRIGATION AND FOOT PEDALS

STERILE IRRIGATION KIT



Sterile irrigation kit.
Box of 5.
Order code: S101461

FOOT PEDALS



- Both wired and wireless foot pedal available.
- Wireless foot pedal operates using a secure encrypted network for complete protection.
- A or B tool selection.
- Foot pedal underlighting to allow for easy detection.

	Model	Order code
Wired	MFP-120	103556
Wireless	MFP-121	103562



4 metre cable for
connecting and
charging foot pedals
via power console.

Order code: 103597

POWER SUPPLY

Mains power cordsets -
3 metres long.
Please specify - Order code

UK	8290
Euro	8300
US/Canada	8310
Australia/New Zealand	8320
Swiss	14700
Brazil	21340



HANDPIECE	Order code	Model	Speed rpm
-----------	------------	-------	-----------



103339 HDS-111 100,000

The HDS-111 handpiece accepts all straight and angled attachments and craniotome

PERFORATOR	Order code	Model	Speed rpm	Operating pressure
------------	------------	-------	-----------	--------------------



103524 P-P02 850 7 bar
Uses Hudson chuck type

FOOT CONTROL	Order code	Model	Foot pedal inlet connection	Foot pedal outlet connection
--------------	------------	-------	-----------------------------	------------------------------



14900 FV-220 Integral 3 metre hose with MA7 probe Accepts 'wash assist' hose from sizes below



3 metre length - Order code 632333
4 metre length - Order code 632343
5 metre length - Order code 632353



Wash cap for use with HDS-111/ P-P02 Order code 9950

STERILISATION TRAYS AND CONTAINER

STERILISATION TRAYS AND CONTAINER



ELECTRIC INSTRUMENT SYSTEM

Order code: 103610

Includes lid, container, electric tray and attachment tray.

Handpiece tray

Optional tray contents -

- 1 off HDS-12X handpiece
- 1 off P-P01 perforator
- 1 off electric cable
- 1 off spanner

Attachment tray

Optional tray contents -

- 5 off attachments
- 1 off spanner

PNEUMATIC INSTRUMENT SYSTEM

Order code 103629

Includes lid, container, pneumatic tray and attachment tray.

Handpiece tray

Optional tray contents -

- 1 off HDS-111 handpiece
- 1 off P-P02 perforator
- 1 off pneumatic hose
- 1 off spanner

Attachment tray

Optional tray contents -

- 5 off attachments
- 1 off spanner

CONTAINER AND TRAYS ONLY

Wire container and lid

Order code: 103944

Tray for electric system

Order code: 103953

Tray for pneumatic system

Order code: 103989

Tray for attachments

Order code: 103966

UNITED KINGDOM (INTERNATIONAL HEAD OFFICE)

De Soutter Medical Limited
Halton Brook Business Park, Weston Road,
Aston Clinton, Aylesbury, Bucks HP22 5WF
Tel: +44 (0) 1296 634 000
Email: info@de-soutter.com

AUSTRALIA

De Soutter Medical Australia Pty Ltd
2/12-14 Apollo Drive, Hallam
Victoria 3803
Tel: +61 (0) 3 9702 4441
Email: australia@de-soutter.com

AUSTRIA

De Soutter Medical Austria
Hauptstrasse 122-124/Büro 01
3420 Kritzensdorf
Tel: +43 (0) 2243 239000
Email: austria@de-soutter.com

BELGIUM

De Soutter Medical Belgium
Bessemersstraat 14
3620 Lanaken
Tel: +32 (0) 89/47 15 37
Email: belgium@de-soutter.com

FRANCE

De Soutter Médical France
1 252 Avenue du Parc des Expositions
33260 La Teste de Buch
Tél. : +33 (0) 5 56 54 89 36
Email : france@de-soutter.com

GERMANY

De Soutter Medical Deutschland
Bahnhofstraße 4
D-66625 Nohfelden
Tel: +49 (0) 68 52-99 12 46
Email: deutschland@de-soutter.com

ITALY

De Soutter Medical Italia
Località Fornace SNC
27022 Casorate Primo - PV
Tel: +39 (0) 2 9009 4098
Email: italy@de-soutter.com

NETHERLANDS

De Soutter Medical Nederland
Gelderlandhaven 2X
3433 PG Nieuwegein
Tel: +31 (0) 85 0491480
Email: nederland@de-soutter.com

UNITED STATES OF AMERICA

De Soutter Medical USA, Inc.
224 Rolling Hill Road, Suite 12A
Mooresville, NC 28117
Tel: +1 (704) 655 9040
Email: usa@de-soutter.com

deSoutter
MEDICAL

Distributed by De Soutter Medical, manufactured by Zethon