

Please visit
our website for
more information




legholderrx

Wirz Medical

Bernstrasse 1
CH-3076 Worb

T +41 31 832 43 53
M +41 79 630 66 44
mechanik@wirz-ortho.ch

legholderrx.ch

CE® mod. dép.

Swissmedic N°: CH201505-0035
GLN Company: 7601 0000 8560 3
Compliant to MDR 13/485
EU Representative: endocon.de
Wirz Ortho Mechanics GmbH
Company-UID (EORI): CHE-480.499.766


legholderrx

A must for knee operations

The **robust foot support** for
6 intraoperative, stable knee positions





Self-retaining knee position with lateral support on the thigh


Facts and technical data

- Easy-care, robust, radiolucent and non-allergenic plastic parts
- Swiss product developed by an orthopedic surgeon and engineer
- Successfully in use in over 40 countries for more than 20 years
- 490x240x200 mm, approx. 2.7 kg



Non-sterile assembly under the drapes
legholderrx.ch

LegholderRX® models



LegholderRX® fusion

The hemispherical nubs on the sides allow the cloth to slide easily between the brackets and the foot is stabilized on the side. | Article no. LRX2024-03

LegholderRX® standard (predecessor model)

Proven and widely used model which has been successfully on the market worldwide for many years. The „standard“ model differs from the „fusion“ model in that it has more-pointed nubs on the sides and slightly thinner brackets. | Article no. LRX2020-02



LegholderRX® extra stable

Improved foot stabilization thanks to biometrically calculated, ergonomic support surfaces and recesses for the Achilles tendon.

| Article no. LRX2024-04

Fastening straps



The specially designed **fastening straps with Velcro fasteners** enable quick and stable installation.

S: Short straps for the majority of operating tables with separate and independently movable leg sections on the right and left.
Advantage: Simpler and more stable assembly with freedom of movement of the individual leg sections.

| Article no. LRX2020-10

L: Long straps are also suitable for operating tables without separate leg sections.

| Article no. LRX2020-11

